
## Quality Assurance Certificate

This certifies that the designated product was manufactured by Sartorius Stedim Biotech in accordance with the current Good Manufacturing Practice standards.

This product has been subjected to and has fulfilled Sartorius Stedim Biotech' rigorous quality control standards from the raw material to the final product.

The Quality Management Department for quality control of membranes, Sartorius Stedim Biotech GmbH, has measured | controlled the values and assures that they are within the limits that are established in the current specification for this material.

Each manufacturing lot was sampled, tested and released by Quality Assurance with respect to the following characteristics,

- Thickness
- Flow Rate
- Bubblepoint visuell
- Growth (pore size < 0.65 µm)

Details of the methodologies used can be obtained from Sartorius Stedim Biotech GmbH.

	Flow rate for water (ml/min/cm <sup>2</sup> /bar)	Bubble point (bar) acc. DIN 58355	Thickness (µm)	Bacteria retention Concentration   filter area (10/7/cm <sup>2</sup> )
<b>Cellulosenitrat (113xx)</b>				
0.22 µm	25	4.2	130	LRV > 7
0.45 µm	70	2.4	130	
0.65 µm	130	2	130	
0.8 µm	200	1.4	130	
1.2 µm	200	1	130	
3 µm	430	0.5	130	
5 µm	570	0.5	130	
8 µm	750	0.3	130	

## Cellulose Nitrate Filters

### Directions for Use

Nitrocellulose membrane filters	
<b>Hazard statements</b>	<b>Precautionary statements</b>
H228: Flammable Solid	P 210: Keep away from heat/spark/open flames/hot surfaces. – No smoking. P 243: Take precautionary measures against static discharge.
	UN 3270; ADRIMDGIIATA 4.1
DANGER	

Producer: Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
D-37079 Goettingen, Germany  
Phone: +49.551.308.0

#### 1. Storage

Store the filters in a dust-free place and not above room temperature. Avoid exposing them to sunlight and solvent or other chemical vapors. Caution: Cellulose nitrate filters are highly flammable. Their flashpoint is approx. 200°C. Always keep these filters separated by the discs of interleaving paper that are supplied with the filters.

#### 2. Labeling

The label on the front of the box provides specifications on the type, number per box, diameter, pore size and the lot number. In case you need to make any inquiries, always state the lot number in your correspondence.

#### 3. Packaging and Removal

Blunt-tipped forceps (16625) are required to remove a filter from the filter tray or package. Avoid touching the filters with your fingers since skin perspiration or pressure may leave non-wettable spots on them. Always hold large-diameter filters in two places to prevent creasing them. The filters are always separated by yellow discs of interleaving paper which must be discarded after removal of a filter.

#### 4. Sterilization

Sartorius cellulose nitrate membrane filters (113xx) can be sterilized by various methods:

##### 4.1.1 Autoclaving

Temperature range [°C]	Pressure range [bar]	Time [min]	Condition
121	1	30–60	Wet

Remarks:

The membrane filters should be wetted with water and placed in a filter holder before they are autoclaved.

It is also necessary to assure that the membrane filter remains wetted throughout the autoclaving cycle and cooling process so that no post-vacuum will be generated. When a dry filter is autoclaved, the pore structure will constrict to a certain degree. This will result in a reduction in the flow rate. However, constriction does not adversely affect the microbial retention capacity.

#### Warning!

The temperature may not exceed 121°C; otherwise, this could cause the membrane to disintegrate.

#### 4.1.2 Sterilization with Ethylene Oxide

Cellulose nitrate membrane filters can be sterilized with ethylene oxide.

#### 4.1.3 Sterilization by γ-Irradiation

The cellulose nitrate membrane filters can be sterilized by γ-irradiation (e.g. with 25 kGy).

#### 4.1.4 Sterilization by Chemical Disinfection

For example, place in 3% aqueous formaldehyde for 25 hours.

#### 4.1.5 Ready-to-Use Filters (already sterilized)

The sterilization method applied for the respective membrane products has been validated.

#### 5. Test Method for Sterile Filtration

##### Bubble Point

The bubble point method is used to detect the largest pores in a membrane filter. Wet the membrane and place in a filter holder. Apply pressure to the membrane. The bubble point is reached when air passes through the largest pores (constant stream of air bubbles).

Remarks:

A certain quantity of air will pass through the membrane before the bubble point is reached. This passage of air is caused by diffusion. As the pressure increases, so will the rate of diffusion. The bubble point depends on the medium used to wet the membrane.

#### 6. Chemical Compatibility

Please refer to our Sartorius Laboratory Product Catalogue.

#### 7. Special Information

Make sure to keep the membrane filter in the same position as in the tray when placing the filter in a filter holder. In other words, the visible surface (upper side as in the tray) must face the inlet of the filter holder in order to obtain the full benefits of the filter.

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Specifications subject to change without notice. Printed and copyrighted by Sartorius Stedim Biotech GmbH