

LABORHITS

Microliter Pipettes LLG-proMLP

Precise and Reliable Pipette

Extremely light weight while the robustness still remains


- High accuracy, reliability and reproducibility
- Perfect ergonomics for highest comfort
- Colour-coded push-button for quick operation
- Extremely low pipetting forces
- Built-in volume locking for error-free results
- Tip ejector button*, adjustable for right- and left-handed users
- Completely autoclavable at 121 °C
- Including individual calibration certificate (according to EN ISO 8655)



*not available for 5 ml and 10 ml pipettes

Capacity µl	Increments µl	Accuracy max. vol. (≤ ± R%)	Precision max. vol. (≤ VK%)	Colour code	PK	Cat No.
0,1 ... 2	0,002	1,65	0,75	dark orange	1	6.270 420
0,5 ... 10	0,02	1,10	0,44	red	1	6.270 421
2 ... 20	0,02	1,00	0,33	light yellow	1	6.270 422
10 ... 100	0,2	0,80	0,17	light orange	1	6.270 423
20 ... 200	0,2	0,80	0,17	dark yellow	1	6.270 424
100 ... 1000	2	0,80	0,17	dark blue	1	6.270 425
500 ... 5000	2	0,66	0,18	purple	1	6.270 426
100 ... 10000	20	0,60	0,18	light blue	1	6.270 427

NEW



LLG-Filter tips Low Retention 2.0, PP

Available in easy-to-open boxes with color coded racks. Compatible with all common pipette brands. Improved optical clarity and pipette fit. Colored coded rack inserts for easy identification. Graduated, sterile, DNase and RNase free

Capacity µl	Package contents	PK	Cat. No.
0,1 ... 10	10 x 96	960	4.668 784
0,1 ... 10*	10 x 96	960	4.668 785
1 ... 20	10 x 96	960	4.668 786
1 ... 100	10 x 96	960	4.668 787
1 ... 200	10 x 96	960	4.668 788
100 ... 1000	8 x 96	768	4.668 789

* Tips XL

LABORHITS

LLG-Pipette Tip Refill System *ULTRALOW*

- DNase, RNase and pyrogen free
- *Ultralow* retention tips with easy visible volume graduations
- Compatible with common pipette brands
- Easy-to-use lift-and-load refill system
- Refill inserts printed with volume and LOT No.
- Racks are stackable, autoclavable and compatible with automated pipette systems

All *ULTRALOW* LLG-Pipette tips feature an extremely low liquid retention compared to standard tips, due to their hydrophobic inner side. This enables accurate and more precise dispensing of fluids, even viscous ones!

LLG refills are easy to use and allow quick filling of empty racks.

LLG-Pipette tip Refill system *ULTRALOW* is compatible with all major pipette brands, such as LLG-proMLP, Gilson, Eppendorf, Brand, Sartorius, Socorex and others.



LLG-Pipette Tips *ULTRALOW*, Refill System

Please order empty racks separately.

Capacity µl	Package contents	PK	Cat. No.
10	10 x 96	960	6.283 565
200	10 x 96	960	6.283 568
1250	8 x 96	768	6.283 570



Racks for LLG-Pipette Tips *ULTRALOW*, empty

Description	PK	Cat. No.
for 10 µl XL and 200 µl LLG- Pipette Tips <i>ULTRALOW</i> Refill System	5	6.283 566
for 1250 µl LLG- Pipette Tips <i>ULTRALOW</i> Refill System	5	6.283 571



LLG-Pipette Tips *ULTRALOW*, Rack

Capacity µl	Package contents	PK	Cat. No.
10	10 racks of 96 tips	960	6.283 567
200	10 racks of 96 tips	960	6.283 569
1250	8 racks of 96 tips	768	6.283 572