

LABORHITS

LLG-Pipette Tip Refill System *ULTRALOW*

Ultralow Retention

**Environmentally Friendly
Refill System**

**No Transfer
Device Needed**



- DNase, RNase and pyrogen free
- *Ultralow* retention tips with easy visible volume graduations
- Universal fit tips, compatible with common pipette brands
- Easy-to-use lift-and-load refill system
- Refill inserts are imprinted with volume and LOT No. for easy identification
- Racks are stackable, autoclavable and compatible with automated pipette systems

All *ULTRALOW* LLG-Pipette tips feature an extremely low liquid retention compared to standard tips, due to their hydrophobic inner side. This enables accurate and more precise dispensing of fluids, even viscous ones!

LLG refills are easy to use and allow quick filling of empty racks.

LLG-Pipette tip Refill system *ULTRALOW* is compatible with all major pipette brands, such as LLG-proMLP, Gilson, Eppendorf, Brand, Sartorius, Socorex and others.

LABORHITS



LLG-Pipette Tips *ULTRALOW*, Refill System

Please order empty racks separately.

Capacity µl	Package contents	PK	Cat. No.
10	10 x 96	960	6.283 565
200	10 x 96	960	6.283 568
1250	8 x 96	768	6.283 570



Racks for LLG-Pipette Tips *ULTRALOW*, empty

Description	PK	Cat. No.
for 10 µl XL and 200 µl LLG- Pipette Tips ULTRALOW Refill System	5	6.283 566
for 1250 µl LLG- Pipette Tips ULTRALOW Refill System	5	6.283 571



LLG-Pipette Tips *ULTRALOW*, Rack

Capacity µl	Package contents	PK	Cat. No.
10	10 racks of 96 tips	960	6.283 567
200	10 racks of 96 tips	960	6.283 569
1250	8 racks of 96 tips	768	6.283 572