

# LABORHITS

## LLG-Buffer solution pH

Accuracy  $\pm 0,01$  pH

With certificate of analysis

Colour-coded

### Advantages:

- Accuracy  $\pm 0.01$  pH, tested and certified at 25 °C
- Temperature dependence data table printed on the label
- Traceable to NIST
- Shelf-life: 24 months (pH 7 and 10) or 36 months (pH 4) from date of production



### Ready-to-use, for easy and exact calibration of your pH-meters.

The ready-to-use colour coded pH-buffer solutions offer excellent quality, precision and are designed to suit all end user requirements. Supplied in 250 ml or 1 L bottles and conveniently colour coded to eliminate any risk of confusion between the 3 pH-values. Certificates of analysis under [www.llg-labware.com](http://www.llg-labware.com).

**Twin-neck dispensing bottles** with integrated calibration vessel, allows the particularly economical dispensing of solution and prevents contamination.

pH-value at 25 °C	Capacity ml	Colour	Cat. No.
4.00	250	red	9.040 868
7.00	250	yellow	9.040 869
10.00	250	blue	9.040 871

### Bottle, 1 liter

pH-value at 25 °C	Capacity ml	Colour	Cat. No.
4.00	1000	red	9.040 865
7.00	1000	yellow	9.040 866
10.00	1000	blue	9.040 867