


LABORHITS

LLG timer

Built in white board



*Up to 4
measurements
simultaneously!*

**Perfect for high noise
locations**



*Vibration,
sound and
light alarm!*

Find out more to 1 and 2

- Count direction: up/down
- Key press confirmation feature gives an audible indication of data entry
- Requires two AAA batteries (included)
- Individually tested, and calibrated using NIST
- Includes Certificate of Calibration with unique report number

Specifications	LLG Timer with White Board 1	LLG-Timer with Triple Alarms 2
	white board and dry erase marker with eraser tip	
Alarm	audible alarm for each channel	Vibration, light and loud buzzer sound alarms can function at one time, or be independently turned off
Kanal	4 channel, adjustable	1 channel
Display	hours / minutes or minutes / seconds	hours / minutes or minutes / seconds
Attachment	flip-open stand for a vertical view and use the magnetic back	flip-open stand for a vertical view and use the magnetic back
Timer accuracy	±5 seconds/day	±5 seconds/day
Switchable timing capacity	23 hours: 59 minutes / 59 seconds	19 hours: 59 minutes / 99 minutes: 59 seconds
Dimensions (W x D x H):	114.3 x 25.4 x 88.9 mm	45 x 75 x 20 mm (without clip portion)
Cat. No.	6.266 706	6.266 707

LABORHITS

LLG timer

Short period timer

Easy-to-use



Specifications	3 LLG-Short period timer with alarm	Interval timer with alarm 4
Running time	timer and stopwatch from/up to 99 min. and 59 sec.	max. 60 mins.
Alarm	with large display and alarm function	interval timer with alarm
Attachment	with magnet for easy fixing	suitable for standing, hanging and in addition on the back of two magnets
Weight	battery LR 44 inclusive approx. 40 g	metal movement 118 g
Dimensions	(Ø x H) 54 x 16.5 mm	(W x D x H) 65 x 65 x 35 mm
Cat. No.	9.260 150	9.260 020