

LABORHITS

LLG-Serological pipettes, type 3

Especially suited for use in clean rooms

Single packed in easy-to-open plastic/plastic packaging

Ideal for cell culture applications



Features:

- Polystyrene (PS)
- Plugged
- Non-pyrogenic
- Easy to read graduation scales facilitate volume measurement and dosing
- Increased volume capacity through extra graduations
- Colour coded for easy volume identification
- Sterile

Volumen ml	Graduation ml	Colour code	Negative graduations up to	PK	Cat. No.
1	0.01	clear	0,3	1000	6.268 099
2	0.01	clear	0,3	1000	6.268 100
5	0.1	blue	2	200	6.268 101
10	0.1	orange	3	500	6.268 102
25	0.2	red	10	200	6.268 103
50	0.5	purple	10	100	6.268 104