

LABORHITS

LLG-uniINCUB 20 and LLG-uniSHAKER 2 Mini Incubator and optional tumbling shaker



NEW



- **Portable, digital, with fan circulation**
- **Perfect for Hygiene Analysis, Dipslides and Contact plates**
- **Perfect for small Laboratories, Universities and Schools**

LLG-Mini-Incubator uniINCUB 20

The portable LLG-uniINCUB 20 digital incubator is ideally suited for haematology and microbiology applications, such as the cultivation of contact plates for microbiological hygiene analysis. Due to its small footprint and economical price, it is also the perfect incubator for educational institutions and small laboratories. The LLG-uniINCUB 20 incubator features a large internal chamber with a capacity of 20 liter, capable of storing flasks and bottles up to 2 liter. In addition, the incubator includes two adjustable/ removable shelves for increased capacity, as well as an electrical socket. The digital temperature control eliminates the need for external thermometers and repetitive „fine tuning“ of an analog control knob. Optimum temperature distribution throughout the cabinet interior by fan circulation. The LLG-mini tumbling shaker uniSHAKER 2 which can be positioned and used inside the incubator, please order separately.

LLG-Tumbling Shaker uniSHAKER 2

Multi-functional shaker with dimpled mat, perfect for foam-free mixing of blood in centrifuge tubes. Suitable for a variety of tubes and dishes.

Specifications LLG-uniINCUB 20:

Temperature range:	Ambient +1 °C to 60 °C
Accuracy:	0.5 °C
Temperature uniformity:	±1.5 °C
Capacity:	20 l
External dimensions (W x D x H):	335 x 370 x 475 mm
Internal dimensions (W x D x H):	260 x 235 x 325 mm
Weight:	6.5 kg
Power supply:	100-240 V, 50/60 Hz
Warranty:	3 Years

Specifications LLG-uniSHAKER 2:

Tumbling angle:	20°
Speed:	Fixed, 20 min ⁻¹
Max. Load:	1.6 kg
Platform Dimensions (W x D):	152,5 x 202,5 mm
Warranty:	3 Years

Type	Cat. No.
LLG-uniINCUB 20 with EU plug	6.263 570
LLG-uniINCUB 20 with UK plug	6.263 571
LLG-uniSHAKER 2 with EU plug	6.263 640
LLG-uniSHAKER 2 with UK plug	6.263 641