

LABORHITS

LLG Syringe Filter



Sterile or non sterile

- Housing in Polypropylene or Acryl
- Multifunctional Syringe filters equipped with Female LuerLock/Male Luer Lock or Male Luer Slip connections for different applications
- Sterile products single packed in individual hard blister
- Superior pressure stability
- Colour coding and printed description for easy identification

Especially designed to provide efficient filtration of a broad variety of fluids, solvents, aqueous or inorganic solutions.

Our syringe filters cover most applications in HPLC, pharmaceutical, environmental, biotechnology and food and beverage testing laboratories.



Sterile

Membrane	Pore size µm	Diam. mm	Colour	PK	Cat. No.
CA	0.20	13	blue	50	9.055 510
CA	0.20	25	blue	50	9.055 511
CA	0.45	13	yellow	50	9.055 512
CA	0.45	25	yellow	50	9.055 513

Non sterile

Membrane	Pore size µm	Housing	Diam. mm	PK	Cat. No.
CA	0.20	Acrylic, blue	13	500	9.055 500
CA	0.20	Acrylic, blue	25	500	9.055 501
CA	0.45	Acrylic, yellow	13	500	9.055 502
CA	0.45	Acrylic, yellow	25	500	9.055 503
CA	0.80	Acrylic, green	13	500	9.055 504
CA	0.80	Acrylic, green	25	500	7.970 389
Nylon	0.20	PP	13	500	9.055 520
Nylon	0.20	PP	25	500	9.055 521
Nylon	0.45	PP	13	500	9.055 522
Nylon	0.45	PP	25	500	9.055 523
RC	0.20	PP	13	500	9.055 530
RC	0.20	PP	25	500	9.055 531
RC	0.45	PP	13	500	9.055 532
RC	0.45	PP	25	500	9.055 533
PE	0.20	PP	13	500	9.055 540
PE	0.20	PP	25	500	9.055 541
PE	0.50	PP	13	500	9.055 542
PE	0.50	PP	25	500	9.055 543
PTFE	0.20	PP	13	500	7.970 402
PTFE	0.20	PP	25	500	9.055 535
PTFE	0.45	PP	13	500	7.970 385
PTFE	0.45	PP	25	500	6.255 331
PES	0.20	PP	13	500	9.055 524
PES	0.20	PP	25	500	9.055 525
PES	0.45	PP	13	500	9.055 526
PES	0.45	PP	25	500	7.970 213
PVDF	0.20	PP	13	500	7.970 258
PVDF	0.20	PP	25	500	7.970 387
PVDF	0.45	PP	13	500	9.055 534
PVDF	0.45	PP	25	500	7.970 286
GF	0.70	PP	25	500	9.055 550
GF	1.00	PP	25	500	9.055 551
GF	1.20	PP	25	500	9.055 552
GF	3.10	PP	25	500	9.055 553