

# LABORHITS

## LLG-Timer pro und LLG-Dual-Timer



**NEW**



### 2 in 1 LLG-Timer pro

**Real time clock and timer at the same time**

The timer also runs when time is visible.

- Digital stopwatch (counting up)
- Short timer (counting down)
- With memory function
- Stopwatch and countdown up to 23 h 59 min. 59 sec.
- Optical and acoustic alarm when timer expires
- Incl. magnet for mounting, hanging bracket and stand (fold-out)

### 2 in 1 LLG-Dual-Timer

**Two separate timers**

One fixed time can be stored per channel.

- Digital stopwatch (counting up)
- Short timer (counting down)
- With memory function
- Stopwatch with tenths of a second
- Stopwatch and countdown up to 99 h 59 min. 59 sec
- Optical (red LED) and acoustic alarm (2 volume levels) when timer expires
- Incl. magnet for mounting, hanging bracket and stand (fold-out)

#### Specifications

Dimensions (W x H x D):  
Weight:  
Power supply:  
Warranty:

#### LLG-Timer pro

82 x 68 x mm  
82 g  
2 x AAA 1,5 V  
(not included)  
12 months

#### LLG-Dual-Timer

97 x 79 x 21 mm  
95 g  
2 x AAA 1,5 V  
(not included)  
12 months

Type	PK	Cat. No.
LLG-Timer pro	1	<b>6.311 351</b>
LLG-Dual-Timer	1	<b>6.311 350</b>