

LLG-Narrow-Neck Bottles

Very stable design

With graduation

Transparent



PP, Heavy Duty, with screw cap

Ideal for collecting, transporting and storing liquid samples.

Capacity ml	Diam. mm	Thread external diam. mm	Height mm	PK	Cat. No.
4	16.0	13	40	100	6.289 520
8	25.0	19	44	100	6.289 450
30	36.0	19	61	100	6.289 451
60	40.0	19	85	100	6.289 452
125	51.0	23	99	50	6.289 453