

# LABORHITS

## Mini Centrifuges LLG-uniCFUGE 2 and 2/5 & LLG-Centrifuge Tubes 5 ml

Compact and reasonably priced

2 speeds selectable



**LLG-uniCFUGE 2/5**  
with rotor for 5 ml tubes

**LLG-uniCFUGE 2**  
with universal rotor



# LABORHITS

## Mini centrifuges LLG-uniCFUGE 2 and 2/5

### Features

- Speed can be adjusted to 4000 rpm or 6000 rpm
- Powerful motor for fast starting
- Quick stop when the lid is opened

### Scope of delivery:

#### LLG-uniCFUGE 2:

with **Universal rotor** for 6 x 1.5/2.0 ml and 6 x 0.5 ml tubes and 2 times 8 x 0.2 ml PCR-strips.

#### LLG-uniCFUGE 2/5:

with rotor for 4 x 5 ml tubes.

### Specifications

**Max. RCF:** 1200 x g and 2000 x g

**Max. speed:** 4000 rpm and 6000 rpm

#### Max. load LLG-uniCFUGE 2:

6 x 0.5/1.5/2 ml and 2 x 8-tube PCR strips (0.2 ml)

#### Max. load LLG-uniCFUGE 2/5:

4 x 5 ml tubes

**Dimensions (diam. x H):** 140 x 120 mm

**Weight:** 0.7 kg

**Power supply:** 100 - 240 V, 50/60 Hz, 1 A




Description	PK	Cat. No.
LLG-uniCFUGE 2 with universal rotor with EU-plug	1	6.263 510
LLG-uniCFUGE 2 with universal rotor with UK-plug	1	6.263 511
LLG-uniCFUGE 2/5 with rotor for 5 ml Tubes with EU-plug	1	6.263 515
LLG-uniCFUGE 2/5 with rotor for 5 ml Tubes with UK-plug	1	6.263 516

## LLG-Centrifuge tubes, 5 ml, PP

### Features

- Max. 25,000 x g, suitable for high speed centrifugation
- Graduations every 0.25 ml for easy volume estimation
- A flat, frosted cap and large writing pad for sample labeling and identification
- Autoclavable
- Certified free of RNase and DNase



Fits most rotors and racks that hold 15 ml conical tubes, no adapters needed. The LLG-Centrifuge tubes 5 ml, represent the perfect solution for all laboratory applications requiring incubation, centrifugation and sample handling in an intermediate volume range of up to 5 ml.

Description	Package contents	PK	Cat. No.
Clear	1 bag of 250 pieces	250	6.281 113
Colour assortment green, yellow, red, blue	4 bags of 50 pieces (1 colour per bag)	250	6.281 114
Amber	1 Beutel à 250 Stück	250	6.281 115