

# LABORHITS

## Disc Rotator LLG-uniLOOPMIX 2

Very flexible  
clamp fixation system

Easy to clean

Cost-efficient



Also order



**LLG-Centrifuge Tubes ECONOMY**  
(6.270 403, 6.270 404, 6.270 405)



Analogue Disc rotator for gentle and effective sample mixing. Suitable for all biological samples. With easily interchangeable clamps that allow a quick cleaning of the rotator. For a maximum of 48 x 1.5 ml, 16 x 5 or 15 ml, 8 x 50 ml tubes. Also for parallel use of different sizes, can be equipped individually.

**Scope of delivery:**

Disc rotator incl. 12 clamps each for 2 ml, 15 ml and 50 ml tubes

**Specifications**

**Type of movement:** rotating  
**Variable speed:** 0 - 80 rpm  
**Rotation angle of disc:** 360°  
**Permissible shaking weight (incl. attachment):** 2 kg  
**rel. Moisture:** up to max 80 %  
**Working temp.:** +5 bis +40 °C  
**Motor rating input/output:** 20 W  
**Dimensions (B x T x H):** 220 x 300 x 310 mm  
**Power supply:** 100 - 240 V, 50/60 Hz  
**Protection class:** IP21

Type	PK	Cat. No.
LLG-uniLOOPMIX 2 with EU-plug	1	6.263 665
LLG-uniLOOPMIX 2 with UK-plug	1	6.263 666
Clamp for 2 ml centrifuge tube	12	6.263 662
Clamp for 15 ml centrifuge tube	12	6.263 663
Clamp for 50 ml centrifuge tube	12	6.263 664