

LABORHITS

Mini-Ultrasonic Cleaner LLG-uniSONIC 1



**Quick and Efficient Cleaning
without Chemical Additives**

Easy Handling

For convenient cleaning of goggles or smaller laboratory instruments. Eco-friendly, only needs warm water. By ultrasonic sound millions of small air bubbles are generated. Combined with a very high pressure these achieve a fast, efficient and effortless cleaning.

Specifications

Weight: 630 g

Power supply: 220 - 240 V, 50 Hz

Outer Dimensions (W x D x H) approx.: 200 x 115 x 115 mm

Dimensions tank (W x D x H) approx.: 155 x 80 x 40 mm

Features

- Ready to use, no assembly required
- 420 ml stainless steel tank (max. filling volume 300 ml)
- Automatic 3-minute-operation, can be stopped anytime
- Removable lid
- 45 kHz ultrasonic frequency
- 30 Watt cleaning performance

Type	PK	Cat. No.
LLG-uniSONIC 1	1	6.263 560