

LABORHITS

Microliter Pipettes LLG-proMLP

Precise and Reliable Pipette

Extremely light weight while the robustness still remains

Features

- High accuracy, reliability and reproducibility
- Perfect ergonomics for better comfort
- Extremely low pipetting forces
- Built-in volume locking for error-free results
- Completely autoclavable at 121 °C
- Including individual calibration certificate (according to EN ISO 8655)

The new LLG-proMLP single channel microliter pipettes combine high accuracy and reliability with ergonomic and extremely low force pipetting relieving stress from your hand and thumb. This high precision instrument has been built to maintain accuracy for many years making it the ideal choice for all scientific work requiring exact and reproducible pipetting. The LLG-proMLP pipettes are equipped with an adjustable eject button allowing right and left handed users to set the button* according to their preferences.

*not available for 5 ml and 10 ml pipettes



Tip ejector button, adjustable for right- and left-handed users

Capacity µl	Increments µl	Accuracy max. vol. (≤ ± R%)	Precision max. vol. (≤ VK%)	Colour code	PK	Cat No.
0.1-2	0.002	1.65	0.75	dark orange	1	6.270 420
0.5-10	0.02	1.10	0.44	red	1	6.270 421
2-20	0.02	1.00	0.33	light yellow	1	6.270 422
10-100	0.2	0.80	0.17	light orange	1	6.270 423
20-200	0.2	0.80	0.17	dark yellow	1	6.270 424
100-1000	2	0.80	0.17	dark blue	1	6.270 425
500-5000	2	0.66	0.18	purple	1	6.270 426
100-10000	20	0.60	0.18	light blue	1	6.270 427