

# LABORHITS

## LLG-pH Pen, pH tester



### Specifications

**Measuring range:** 0.00 to 14.00 pH

**Accuracy:** ±0.1 pH

**Dimensions with electrode (L x W x H):**

210 x 36 x 20 mm

**Dimensions electrode (diam. x L):**

12 x 110 mm

**Weight:** 72 g

**Battery:** 2 x 3 V CR2032 button cell

**Protection class:** IP65

**Warranty:** 2 years

For measurement in  
narrow vessels

Cost efficient pH tester

3-point calibration



### Also order

#### LLG-pH-Buffer solutions

(9.040 868, 9.040 869, 9.040 871)

### Features

- Ideal for testing small amounts of liquid and in narrow vessels
- 3-point calibration
- ATC from 0 °C and 50 °C
- Supplied with storage vessel, no dry-out of the electrode

Type	PK	Cat. No.
LLG-pH Pen	1	6.263 690