

# LABORHITS

## Universal drying oven LLG-uniOVEN 42 and LLG-uniOVEN 110

42 liter or 110 liter volume

Quick and intuitive  
operation



**NEW**



### Features

- Microprocessor controlled universal drying oven with forced convection
- Digital display for time and temperature
- Timer: 1 minute to 99.9 h and continuous
- Delayed start timer: 1 minute to 99.9 h
- Acoustical and visible information for end of program
- Acoustical and optical alarm at temperature differences
- Outlet for vapour exhaustion
- Interior made of stainless steel
- Adjustable over temperature protection up to 300 °C

**WARRANTY**  
**3**  
**YEARS**

### Specifications

**Internal volume:** 42 l or 110 l  
**Working temperature:** +5 °C up to +40 °C  
**Temperature range:** RT +7 °C up to 250 °C  
**Temperature accuracy:** ± 1 °C (at 70 °C)  
**Temperature uniformity:** ± 2 °C (at 70 °C)  
**Maximum humidity:** 80 %, non-condensing  
**Power supply:** 230 V, 50/60 Hz

**Scope of supply:** LLG-uniOVEN 42 or 110 incl. 2 shelves

Type	Power consumption W	Internal dimensions mm (W x D x H)	External dimensions mm (W x D x H)	Weight kg	PK	Cat. No.
LLG-uniOVEN 42	850	420 x 280 x 360	705 x 505 x 540	39	1	6.263 670
LLG-uniOVEN 110	1550	500 x 445 x 500	790 x 665 x 680	66	1	6.263 675

# LABORHITS

## Universal drying LLG-uniOVEN 150



**LLG-uniOVEN 150**

**WARRANTY**  
**3**  
**YEARS**

**150 liter volume**

**Quick and intuitive operation**



### Specifications

**Internal volume:** 150 l

**Working temperature:** +5 °C up to +35 °C

**Temperature range:** RT +10 °C up to 300 °C

**Temperature accuracy:** ± 0,5 °C (at 70 °C)

**Temperature uniformity:** ± 2 °C (at 70 °C)

**Maximum humidity:** 80 %, non-condensing

**Power supply:** 230 V, 50/60 Hz

### Features

- Microprocessor controlled universal drying oven with forced convection
- Integrated auto-diagnostic system
- 5 different programs with 5 steps can be programmed by the user
- Controller can be protected by password
- Lockable door handle
- LED display
- Interior made of stainless steel (AISI-304)
- Fan speed adjustable
- Hot indicator for temperatures above 50 °C
- Adjustable over temperature protection

**Scope of supply:** LLG-uniOVEN 150 incl. 2 shelves

Type	Power consumption W	Internal dimensions mm (W x D x H)	External dimensions mm (W x D x H)	Weight kg	PK	Cat. No.
<b>LLG-uniOVEN 150, with EU plug</b>	3000	495 x 440 x 705	700 x 665 x 1010	85	1	<b>6.263 680</b>
<b>LLG-uniOVEN 150, with UK plug</b>	3000	495 x 440 x 705	700 x 665 x 1010	85	1	<b>6.263 682</b>