

# LABORHITS

## LLG Laboratory Chairs and Stools

Revolving, height adjustable seat

Easy cleaning and disinfecting


Metal 5-star bar



**1**  
LLG Chair PU  
Standard II plus



**2**  
LLG Chair PU  
Standard II



**3**  
LLG Chair PU  
Standard I plus



**4**  
LLG Stool PU  
Special plus



**5**  
LLG Stool PU  
Standard plus



**6**  
LLG Stool PU  
Standard

- According to EN 1728, EN 1021-1 and 1021-2, EN 61340-5-1.
- Metal 5-star bar with 5 castors for soft and 5 castors for hard and soft floors or gliders.
- Armrests and chromed footring as outlined below.

### LLG Chair PU Foam:

Angle and height of the backrest can be adjusted which always guarantees an ergonomic seating position.

Type	Seat height adjustment mm	Description	dimensions Seat (L x W) mm	PK	Cat. No.
<b>1</b> Chair PU Standard II plus	525 - 770	5 gliders + footring + armrests	440 x 460	1	<b>9.732 200</b>
<b>2</b> Chair PU Standard II	450 - 580	2 x 5 castors, soft floors + 5 castors, hard floors	440 x 460	1	<b>9.732 201</b>
<b>3</b> Chair PU Standard I plus	525 - 770	5 gliders + footring	420 x 400	1	<b>9.732 202</b>
<b>4</b> Stool PU Special plus	525 - 770	5 gliders + footring	350 x 300	1	<b>9.732 204</b>
<b>5</b> Stool PU Standard plus	525 - 770	5 gliders + footring	340 x 340	1	<b>9.732 208</b>
<b>6</b> Stool PU Standard	425 - 540	2 x 5 castors for hard and soft floors	340 x 340	1	<b>9.732 209</b>