

# LABORHITS

## uniPLATESHAKER 1

**Variable speed from  
300 to 1800 rpm**

**Suitable for use  
in cold-rooms**

**Very small  
footprint**



### Features

- Universal fit, comes with interchangeable clamps for different micro plate sizes
- Digital Display with Speed and Time settings
- Maintenance free brushless DC motor
- Orbital Diameter of 2 mm for thorough mixing
- Timer setting from 1 to 99 min and continuous mode
- Comes with foam attachment for 20 x 1.5/2 ml microtubes
- 3-years warranty

The LLG-uniPLATESHAKER 1 has been designed for accommodating and mixing a single micro plate. Due to its small footprint it can be used in working areas with limited space such as fume hoods and bio safety cabinets.

### Specifications

- Type of movement:** orbital
- Orbital diameter:** 2 mm
- Speed range:** 300 - 1800 rpm
- Motor:** Brushless DC motor
- Timer:** 1 - 99 min, continuous
- Power consumption:** 15 W
- Permissible ambient temperature:** 4 - 40 °C
- Dimensions (W x D x H):** 127 x 85 x 86 mm
- Weight:** 1.4 kg
- Power supply:** 100 - 230 V, 50/60 Hz

Type	Cat. No.
uniPLATESHAKER 1 with EU-plug	6.263 488
uniPLATESHAKER 1 with UK-plug	6.263 489