

LABORHITS



**ergonomic
soft grip**

**colour-coded
volume**

LLG - Variable Volume Microlitre Pipettes

Ergonomic shape for a very comfortable handling

The rounded shape and the unique soft grip of the LLG pipette ensures a safe and comfortable handling. A minimal effort is required for pipetting so that the risk of repetitive strain injuries (RSI) is reduced.

Maximum Precision and Accuracy

Thermal insulation of the internal components, using the innovative soft grip, improves accuracy.

Universal Blow-out System

The two-step plunger operation allows reverse pipetting. This also facilitates last drop dispensing by soft plunger movements. High accuracy and precision are guaranteed.

Colour-coded push-button for quick operation

The colour code enables a quick choice of the right pipette depending on the volume to be pipetted (see corresponding column in the table below).

Easy Maintenance and Adjustment

The volume setting is easy and precise, using the push-button, even when wearing gloves. Maintenance and adjustment take only a few minutes with just one tool. (Included in delivery)

Contamination-free Pipetting

The pipette is fully autoclavable, ideal for applications where sterility is important.

Satisfaction Guaranteed

Each pipette is individually tested and calibrated twice, conforming to EN ISO 8655 and DIN 12600 Standards. A certificate is included.

Made in Germany, CE certified.



Available with fixed or adjustable volume.

Capacity µl	Accuracy ≤ ± A% (max. vol.)	Reproducibility ≤ ± CV% (max. vol.)	Colour code	PK	Cat. No.
0.1 to 2,5	1.5	0.8	grey	1	9.280 000
0.5 to 10	1.0	0.5	grey	1	9.280 001
2 to 20	4.5	1.0	grey	1	9.280 002
10 to 100	2.8	0.8	yellow	1	9.280 003
20 to 200	2.5	0.6	yellow	1	9.280 004
100 to 1000	2.9	0.6	blue	1	9.280 005
500 to 5000	2.4	0.6	purple	1	9.280 006
1000 to 10.000*	3.0	0.6	turquoise	1	9.280 007

* 10.000µl Pipettor does not have a tip ejector