

# LABORHITS

## Pipette Tips and Filter Tips „Premium“



metal-free



**Made in Germany**

- **Secure fit and simple ejection of tip**

Optimized cone geometry with wide sealing zone

- **Quick and simple volume check**

Graduation marks on tip

- **Precise pipetting without loss of sample**

Conical tip outlet prevents drop formation at the outlet opening of tip

- **Eliminates the risk of sample contamination**

Tips manufactured from autoclavable, virgin polypropylene

- **Maintain sample integrity**


Available as sterilized version (by beta radiation)

- **Available in bags or rackboxes**



- **No sample contamination from filter**

Use of virgin high quality hydrophobic polyethylene filters without additives




- **No pipette contamination**

Restrains aerosols and liquids (vaporous, radioactive, biohazardous or corrosive samples) through micropore filter technology (MPS technology, filter pore size 8-10 µm)

- **Maximum sample recovery**

**Minimal sample loss**  
**High-precision pipetting**  
Innovative Low Retention surface



### Specifications Filter Tips:

General requirements  
CE-marking certificate









Sterility certificate  
RoHS Compliance certificate  
Environmental Management certificate  
USP-Plastic Test, Class VI, certificate  
Conformity testing certificate  
GMP/GLP certificate

ISO Standard 8655  
EU-Directive 98/79/EG, competent authority code DE/CA90  
EN 552/ISO 11137 and DIN EN 556  
EU Directive 2002/95/EC  
ISO Standard 14000 and ISO Standard 14001  
ISO Standard 10993  
DIN Standard 12600  
WHO Technical Report Series 908, 37. Report, Geneva 2003

# LABORHITS

## Pipette Tips and Filter Tips „Premium“

Images in original size

Pipette tips, low retention	Length	PK	Cat. No.
 0.1-10 µl, clear	34.1 mm	1000 10 x 96	9.409 038 9.409 039
 1-200 µl, clear	50.8 mm	1000 10 x 96	9.409 041 9.409 042
 100-1000 µl, long, clear	86.0 mm	1000 10 x 96	9.409 043 9.409 044
Filter tips, low retention	Length	PK	Cat. No.
 0.1-10 µl, clear, sterile	34.1 mm	10 x 96	9.409 033
 1-20 µl, clear, sterile	50.8 mm	10 x 96	9.409 034
 1-100 µl, clear, sterile	50.8 mm	10 x 96	9.409 035
 1-200 µl, clear, sterile	50.8 mm	10 x 96	9.409 036
 100-1000 µl, long, clear, sterile	86.0 mm	10 x 96	9.409 037