

LABORHITS

LLG-Electrical pipette pump *Plus*

Ergonomic

hand-neutral design
concave padded triggers
lightweight

Very precise and reliable

thanks to microprocessor

Flexible

easy and sensitive control of the aspiration and dispensing speed
5 colored nose pieces



Specifications:

- For 1-100ml glass or plastic pipettes
- Soft touch triggers
- 5 differently colored autoclavable nose pieces
- Uses standard 0.45µm hydrophobic filters
- Long-life Li-Ion battery: up to 20 h continuous pipetting
- Sensitive and adjustable aspiration and dispensing
- Bright LED's to indicate the used mode: High, Low or Gravity Mode
- Labelling area for identification
- Universal voltage and plug faced charger, bench stand and magnetic/fixable wall mount in delivery included



Description	PK	Cat. No.
LLG-Electrical pipette pump <i>Plus</i>	1	9.281 041
Membrane filter, 0.45µm, PTFE, non-sterile	5	9.281 005