

LABORHITS

Mini-centrifuge **LLG-uniCFUGE 5**

Easy switching between RPM and RCF

Large LCD Display

Maximum 15,595 x g

Advantages

- Microprocessor controlled brushless motor for speeds up to 15,000 rpm
- Short-spin mode with selectable RZB and RPM values
- Fixed angle rotor with lid, made of aluminium for 12 x 1.5/2 ml tubes
- Intuitive and userfriendly interface
- Integrated imbalance-detection with automatic shutdown
- Wide suction feet provide a stable base
- The centrifuge can be controlled via computer
- 3 years warranty

Item consists of:

- Centrifuge with rotor (incl. lid) for 12 x 1.5/2 ml tubes
- Power supply with USB cable
- 12 Adapters for 0.4/0.5 ml tubes
- 12 Adapters for 0.2 ml tubes
- Software CD
- User manual

Optional:

- Rotor (Cat. No. 6.263 535) for 16 x 0.2 ml tubes or 2 x 8 PCR-strips



Specifications

Max. RCF:	15,595 x g
Max. speed:	variable from 500 up to 15,000 rpm
Accuracy speed:	±100 rpm
Motor:	brushless DC motor
Max. load:	12 x 2 ml PCR tubes
Dimensions (WxDxH):	262 x 230 x 131 mm
Weight:	4 kg
Power supply:	100 - 240 V, 50/60 Hz

Type	Cat. No.
LLG-uniCFUGE 5 with EU-plug	6.263 530
LLG-uniCFUGE 5 with UK-plug	6.263 531
Rotor for 16 x 0.2 ml tubes or 2 x 8 PCR-strips	6.263 535