

LABORHITS

LLG-uniDISPENS pro

NEW

120 Dosing Programs

Ergonomic

Precise

WARRANTY
3
YEARS



LLG-uniDISPENS pro manual Serial Dispenser

The **ergonomically designed** manual dispenser LLG-uniDISPENS pro enables single-handed volume selection, loading and dispensing. Due to its lightweight design, it is ideally suited for serial dosing.

You can flexibly select your desired dosing volume from **120 dosing programs** for 12 tip sizes and with 10 stroke adjustments. You can view all settings and volumes at a single glance through the viewing window of the manual dispenser. The dispenser is **compatible with Ritips[®] from Ritter, PD-Tips[®] from Brand, Combitips[®] and Combitips advanced[®] from Eppendorf[™].**

Thanks to a purely mechanical, maintenance-free construction, the manual dispenser is a practical and convenient device for daily laboratory use.

Dosing volume 1 µl to 5500 µl.

Scope of supply:

Dispenser, adapter for 25 ml tips

Type	PK	Cat. No.
LLG-uniDISPENS pro, EU-Stecker	1	9.284 239