

LABORHITS


LLG-uniTOPDISPENS

NEW!

Accurate and reliable volume setting

Valve block made of PP gives optimised protection of the glass made dosing cylinder

Precise grinding of the glass piston into the glass cylinder forms a perfect seal



Delivered with 2 PP-adapters
GL 40 and 45



Visible dosing procedure

Ergonomic finger guard

MADE IN GERMANY

LLG-Bottletop dispenser uniTOPDISPENS with glass piston and clear glass cylinder

A reasonably priced bottle-top dispenser with PTFE-coated glass piston available in 2 volumes: 10ml and 50ml. Accurate, reliable and robust. The PTFE-coating of the glass piston is ideal when dispensing of crystallizing media. The high quality of the material used and the precision of the grinding ensures the perfect seal of the complete dispensing system.

Volume ml	Graduations ml	Accuracy max. vol. ($\leq \pm R\%$)	Precision max. vol. ($\leq CV\%$)	PK	Cat. No.
2 - 10	0.5	0.6	0.2	1	9.283 806
10 - 50	1.0	0.6	0.2	1	9.283 808